

MID YEAR SALE

UNBEATABLE PRICES
ON ELECTRONIC AND
MECHANICAL PIPETTES

Hurry, offer ends 29th July 2022.



Call us now to place your order 1800 069 161
or visit [pathtech.com.au](https://www.pathtech.com.au)
for more info or to buy online

 **Pathtech**[™]
First Choice In Science




ELECTRONIC PIPETTES

The Picus Pipette family

REDUCE YOUR RSI RISK

Enjoy your work with the award winning Picus! Its highly sophisticated and completely ergonomic design will reduce your risk of repetitive strain injury (RSI). Fully electronic operation ensures the maximum accuracy and repeatability of your results.

- › Uniquely lightweight and streamlined design for effortless pipetting experience
- › Conveniently located soft-touch operating button and electronic tip ejection help minimize muscle strain
- › The comfortable design of the handle and finger hook enable Picus to rest lightly in your hand, with minimum effort
- › Electronic piston control and brake guarantee accurate and precise pipetting results, independent of the user

 <div style="background-color: red; color: white; border-radius: 50%; width: 40px; height: 40px; display: flex; align-items: center; justify-content: center; margin: 10px auto;"> SAVE \$351 </div>	Picus single channel electronic pipette	Product Code	Original Price	SALE PRICE
	0.2 – 10ul	171-735021	\$1,250.00	\$899.00
	5 – 120ul	171-735041	\$1,250.00	\$899.00
	10 – 300ul	171-735061	\$1,250.00	\$899.00
	50 – 1000ul	171-735081	\$1,250.00	\$899.00
	100 – 5000ul	171-735101	\$1,250.00	\$899.00
500 – 10,000ul	171-735111	\$1,250.00	\$899.00	
 <div style="background-color: red; color: white; border-radius: 50%; width: 40px; height: 40px; display: flex; align-items: center; justify-content: center; margin: 10px auto;"> SAVE \$591 </div>	Picus 8 channel electronic pipette	Product Code	Original Price	SALE PRICE
	0.2 - 10ul	171-735321	\$2,070.00	\$1,479.00
	5 - 120ul	171-735341	\$2,070.00	\$1,479.00
	10 - 300ul	171-735361	\$2,070.00	\$1,479.00
	50 - 1200ul	171-735391	\$2,070.00	\$1,479.00
 <div style="background-color: red; color: white; border-radius: 50%; width: 40px; height: 40px; display: flex; align-items: center; justify-content: center; margin: 10px auto;"> SAVE \$691 </div>	Picus 12 channel electronic pipette	Product Code	Original Price	SALE PRICE
	0.2 - 10ul	171-735421	\$2,390.00	\$1,699.00
	5 - 120ul	171-735441	\$2,390.00	\$1,699.00
	10 - 300ul	171-735461	\$2,390.00	\$1,699.00
	50 - 1200ul	171-735491	\$2,390.00	\$1,699.00

Terms & Conditions: All prices exclude GST. Sale Prices valid on listed products only. Offer ends 29th July 2022.




ELECTRONIC PIPETTES

The Picus® NxT Pipette family

THE PREMIUM ELECTRONIC PIPETTE

The Picus® NxT Electronic Pipette brings ergonomics and advanced function to a new level for your pipetting needs. Picus® NxT takes on where Picus® left with added features like user customizable protocols, password protection, and calibration reminders. Especially when conformity to the strictest regulations is needed by selecting the advanced Picus® NxT you can take your compliance to a new level with added certification and safety features. All the benefits of the Picus plus:

- › Highly advanced, compact and lightweight
- › Advanced applications, modes and protocols
- › Certificate of accredited 3-point calibration to meet demanding quality system requirements
- › RFID-enabled to support asset management

 SAVE \$292	Picus NXT single channel pipette	Product Code	Original Price	SALE PRICE
	0.2-10ul	171-LH-745021	\$1,491.00	\$1,199.00
	5-120ul	171-LH-745041	\$1,491.00	\$1,199.00
	10-300ul	171-LH-745061	\$1,491.00	\$1,199.00
	50-1000ul	171-LH-745081	\$1,491.00	\$1,199.00
	100-5000ul	171-LH-745101	\$1,491.00	\$1,199.00
500-10000ul	171-LH-745111	\$1,491.00	\$1,199.00	
 SAVE \$461	Picus NXT 8 channel pipette	Product Code	Original Price	SALE PRICE
	0.2-10ul	171-LH-745321	\$2,460.00	\$1,999.00
	5-120ul	171-LH-745341	\$2,460.00	\$1,999.00
	10-300ul	171-LH-745361	\$2,460.00	\$1,999.00
	50-1200ul	171-LH-745391	\$2,460.00	\$1,999.00
 SAVE \$538	Picus NXT 12 channel pipette	Product Code	Original Price	SALE PRICE
	0.2-10ul	171-LH-745421	\$2,837.00	\$2,299.00
	5-120ul	171-LH-745441	\$2,837.00	\$2,299.00
	10-300ul	171-LH-745461	\$2,837.00	\$2,299.00
	50-1200ul	171-LH-745491	\$2,837.00	\$2,299.00

Terms & Conditions: All prices exclude GST. Sale Prices valid on listed products only. Offer ends 29th July 2022.

Order today or call 1800 069 161 for more information.




MECHANICAL PIPETTES

The Mline Pipette family

WHERE PRECISION MEETS ACCURACY

Sartorius' mechanical pipette family - the Mline - offers excellent ergonomics, performance and safety in manual pipetting. It is particularly designed to maintain high accuracy and precision in repetitive and long-lasting pipetting. The Mline is an advanced and ergonomic choice for quality-conscious users who perform long and repetitive manual pipetting and is recommended by health and safety officers all over the world.

- ▶ Excellent ergonomics to protect from injuries with the exceptionally light pipetting and tip ejection forces
- ▶ High accuracy and precision, also in repetitive, long-lasting pipetting
- ▶ Minimised risk of contamination with Safe-Cone Filters and full autoclavability
- ▶ Increased safety with volume lock preventing accidental volume changes while pipetting

 SAVE \$96	Mline single channel pipette 0.1 – 3ul 0.5 – 10ul 2 – 20ul 10 – 100ul 20 – 200ul 100 – 1000ul 500 – 5000ul 1– 10ml	Product Code 171-725010 171-725020 171-725030 171-725050 171-725060 171-725070 171-725080 171-725090	Original Price \$375.00 \$375.00 \$375.00 \$375.00 \$375.00 \$375.00 \$375.00 \$375.00	SALE PRICE \$279.00 \$279.00 \$279.00 \$279.00 \$279.00 \$279.00 \$279.00 \$279.00
 SAVE \$166	Mline 8 channel pipette 0.5 - 10ul 5 - 100ul 30 - 300ul	Product Code 171-725120 171-725130 171-725140	Original Price \$795.00 \$795.00 \$795.00	SALE PRICE \$629.00 \$629.00 \$629.00
 SAVE \$206	Mline 12 channel pipette 0.5 - 10ul 5 - 100ul 30 - 300ul	Product Code 171-725220 171-725230 171-725240	Original Price \$975.00 \$975.00 \$975.00	SALE PRICE \$769.00 \$769.00 \$769.00

Terms & Conditions: All prices exclude GST. Sale Prices valid on listed products only. Offer ends 29th July 2022.




MECHANICAL PIPETTES

The Tacta Pipette family

THE ULTIMATE IN ERGONOMICS

Tacta has been designed to feel great in your hands. Thanks to its ergonomic handle and finger hook, there is no need to grip the pipette tightly. Easy, safe and comfortable to use, Tacta meets the most demanding pipetting requirements, delivering consistently accurate results time after time.

- › Comfortable to hold due to the ergonomically designed handle – Weight only: 70-100g single-channel models
- › Low pipetting and tip ejection forces that reduce the risk of RSI and Work Related Upper Limb Disorder
- › The unique Sartorius Optilock system provides flexibility for volume adjustment and locking
- › Controlled and smooth tip ejection with the new Sartorius Optiject technology (levered tip ejection system)

 SAVE \$61	Tacta single channel pipette 0.1 – 3ul 0.5 – 10ul 2 – 20ul 10 – 100ul 20 – 200ul 100 – 1000ul 500 - 5000ul 1– 10ml	Product Code 171-LH-729010 171-LH-729020 171-LH-729030 171-LH-729050 171-LH-729060 171-LH-729070 171-LH-729080 171-LH-729090	Original Price \$410.00 \$410.00 \$410.00 \$410.00 \$410.00 \$410.00 \$410.00 \$410.00	SALE PRICE \$349.00 \$349.00 \$349.00 \$349.00 \$349.00 \$349.00 \$349.00
 SAVE \$236	Tacta 8 channel pipette 0.5 - 10ul 5 - 100ul 30 - 300ul	Product Code 171-LH-729120 171-LH-729130 171-LH-729140	Original Price \$1,085.00 \$1,085.00 \$1,085.00	SALE PRICE \$849.00 \$849.00 \$849.00
 SAVE \$281	Tacta 12 channel pipette 0.5 - 10ul 5 - 100ul 30 - 300ul	Product Code 171-LH-729220 171-LH-729230 171-LH-729240	Original Price \$1,320.00 \$1,320.00 \$1,320.00	SALE PRICE \$1,039.00 \$1,039.00 \$1,039.00

Terms & Conditions: All prices exclude GST. Sale Prices valid on listed products only. Offer ends 29th July 2022.

Order today or call 1800 069 161 for more information.




MECHANICAL PIPETTES

The Proline Plus Pipette family

THE PREMIUM STUDENT PIPETTE

Proline Plus mechanical pipette family is designed to offer comfort and quality for your everyday manual pipetting. It combines durable construction with ease and lightness of use and is therefore a perfect choice for liquid-handling professionals and students alike. Proline Plus has low pipetting forces, a comfortable handle and a robust feel for heavier use.

- ▶ The widest range of pipettes for various users and applications
- ▶ Ergonomic design with light pipetting forces, comfortable handle and finger support for reduced risk of strain injuries
- ▶ Minimised risk of contamination with Safe-Cone Filters and full autoclavability
- ▶ Highly durable with strengthened structure also for heavier use

 SAVE \$76	Proline Plus single channel pipette 0.1-3ul 0.5-10ul 2-20ul 5-50ul 10-100ul 20-200ul 100-1000ul 500-5000ul 1-10ml	Product Code 171-728010 171-728020 171-728030 171-728040 171-728050 171-728060 171-728070 171-728080 171-728090	Original Price \$295.00 \$295.00 \$295.00 \$295.00 \$295.00 \$295.00 \$295.00 \$295.00 \$295.00	SALE PRICE \$219.00 \$219.00 \$219.00 \$219.00 \$219.00 \$219.00 \$219.00 \$219.00 \$219.00
 SAVE \$211	Proline Plus 8 channel pipette 0.5-10ul 10-100ul 30-300ul	Product Code 171-728120 171-728130 171-728140	Original Price \$760.00 \$760.00 \$760.00	SALE PRICE \$549.00 \$549.00 \$549.00
 SAVE \$336	Proline Plus 12 channel pipette 0.5-10ul 10-100ul 30-300ul	Product Code 171-728220 171-728230 171-728240	Original Price \$975.00 \$975.00 \$975.00	SALE PRICE \$639.00 \$639.00 \$639.00

Terms & Conditions: All prices exclude GST. Sale Prices valid on listed products only. Offer ends 29th July 2022.

MECHANICAL PIPETTES

The Proline Pipette family

IDEAL FOR TEACHING LABS

The fact that Proline, Sartorius' first mechanical pipette, is still in use in many laboratories by over 200,000 users, testifies to its timeless, practical design and reliability. Being the most affordable pipette in Sartorius' range of mechanical pipettes, it is ideal for universities and colleges, or any laboratory seeking a cost-efficient liquid handling tool.

- › Wide range of adjustable pipette volumes
- from 2.5ul to 5 ml
- › Easy volume setting with click stop mechanism
- › Safe-Cone Filters available for models > 10ul
- › Relatively lightweight, highly accurate and precise

	Proline single channel pipette	Product Code	Original Price	SALE PRICE
	0.1 – 2.5ul	171-720010	\$180.00	\$139.00
	0.5 – 10ul	171-720015	\$180.00	\$139.00
	2 – 20ul	171-720080	\$180.00	\$139.00
	5 – 50ul	171-720025	\$180.00	\$139.00
	10 – 100ul	171-720050	\$180.00	\$139.00
	20 – 200ul	171-720070	\$180.00	\$139.00
	100 – 1000ul	171-720060	\$180.00	\$139.00
	1 – 5ml	171-720110	\$180.00	\$139.00

Terms & Conditions: All prices exclude GST. Sale Prices valid on listed products only. Offer ends 29th July 2022.



Pipette Service & Calibration

Pathtech provides calibration, maintenance and repair services to extend the lifetime and ensure the continuous accuracy and precision of your pipettes via our NATA Accredited Service Laboratory.

Why is Calibration and Maintenance Needed?

All pipette manufacturers recommend that regular maintenance and calibration is applied to maintain reliable pipetting results and to maximize the lifetime of the pipette. Pipette calibration is a fundamental part of Good Laboratory Practice (GLP) and quality systems and must be considered a vital part of any laboratory regime where precise volumes of fluid need transferring or diluting. Pipette performance is measured as accuracy and precision or how close the dispensed volumes are to the target and how close the results are to each other.

Contact our Pipette Service Team today on 1800 069 161 or email info@pathtech.com.au

Call us now to place your order 1800 069 161
or visit pathtech.com.au for more info or to buy online

Call us now to place your order 1800 069 161
or visit pathtech.com.au
for more info or to buy online



Pathtech[™]
First Choice In Science