

Mid-Year Pipette Promotion

ON SALE UNTIL 31ST JULY 2024



**UNBEATABLE PRICES ON ELECTRONIC
AND MECHANICAL PIPETTES**

Call us now to place your order 1800 069 161
or visit pathtech.com.au for more info or to buy online

SARTORIUS

 **Pathtech**[™]
First Choice In Science

ELECTRONIC PIPETTES


The Picus 2 Pipette family


CONNECTED FOR THE FUTURE


Picus 2 pipettes ensure reliable, repeatable pipetting results and feature an unbeatable ergonomic design that is kind to your hand. Picus 2 is suited to all lab workers, from graduate students to experienced laboratory professionals. It is as intuitive to use as a mechanical pipette but offers more advanced options for experienced users.

By connecting Picus 2 to your mobile device, you can use the Sartorius Pipetting mobile app to smoothly run sample preparation workflows and adjust the pipette setting automatically, taking your productivity to the next level. The app also lets you update the pipettes to the latest software.

Experience the Picus 2 by using the 360-degree viewer tool below to further understand an in-depth overview of features and capabilities.

	Picus 2 - single channel electronic pipette	Product Code	Original Price	SALE PRICE
SAVE \$438	0.2-10ul	171-LH-747021	\$1,437.00	\$999.00
	5-120ul	171-LH-747041	\$1,437.00	\$999.00
	10-300ul	171-LH-747061	\$1,437.00	\$999.00
	50-1000ul	171-LH-747081	\$1,437.00	\$999.00
	100-5000ul	171-LH-747101	\$1,437.00	\$999.00
	500-10,000ul	171-LH-747111	\$1,437.00	\$999.00

	Picus 2 - 8 channel electronic pipette	Product Code	Original Price	SALE PRICE
SAVE \$736	0.2-10ul	171-LH-747321	\$2,435.00	\$1,699.00
	5-120ul	171-LH-747341	\$2,435.00	\$1,699.00
	10-300ul	171-LH-747361	\$2,435.00	\$1,699.00
	50-1200ul	171-LH-747391	\$2,435.00	\$1,699.00

	Picus 2 - 12 channel electronic pipette	Product Code	Original Price	SALE PRICE
SAVE \$772	0.2-10ul	171-LH-747421	\$2,771.00	\$1,999.00
	5-120ul	171-LH-747441	\$2,771.00	\$1,999.00
	10-300ul	171-LH-747461	\$2,771.00	\$1,999.00
	50-1200ul	171-LH-747491	\$2,771.00	\$1,999.00

Terms & Conditions: All prices exclude GST. Sale Prices valid on listed products only. Offer ends 31st July 2024.

SARTORIUS

 **Pathtech**[™]
First Choice In Science


MECHANICAL PIPETTES


The Tacta Pipette family


THE ULTIMATE IN ERGONOMICS

Tacta has been designed to feel great in your hands. Thanks to its ergonomic handle and finger hook, there is no need to grip the pipette tightly. Easy, safe and comfortable to use, Tacta meets the most demanding pipetting requirements, delivering consistently accurate results time after time.

- ▶ Comfortable to hold due to the ergonomically designed handle
 - Weight only: 70-100g single-channel models
- ▶ Low pipetting and tip ejection forces that reduce the risk of RSI and Work Related Upper Limb Disorder
- ▶ The unique Sartorius Optilock system provides flexibility for volume adjustment and locking
- ▶ Controlled and smooth tip ejection with the new Sartorius Optiject technology (levered tip ejection system)

	Tacta single channel pipette	Product Code	Original Price	SALE PRICE
SAVE \$125	0.1 - 3ul	171-LH-729010	\$524.00	\$399.00
	0.5 - 10ul	171-LH-729020	\$524.00	\$399.00
	2 - 20ul	171-LH-729030	\$524.00	\$399.00
	10 - 100ul	171-LH-729050	\$524.00	\$399.00
	20 - 200ul	171-LH-729060	\$524.00	\$399.00
	100 - 1000ul	171-LH-729070	\$524.00	\$399.00
	500 - 5000ul	171-LH-729080	\$524.00	\$399.00
	1 - 10ml	171-LH-729090	\$524.00	\$399.00

	Tacta 8 channel pipette	Product Code	Original Price	SALE PRICE
SAVE \$312	0.5 - 10ul	171-LH-729120	\$1,311.00	\$999.00
	5 - 100ul	171-LH-729130	\$1,311.00	\$999.00
	30 - 300ul	171-LH-729140	\$1,311.00	\$999.00

	Tacta 12 channel pipette	Product Code	Original Price	SALE PRICE
SAVE \$356	0.5 - 10ul	171-LH-729220	\$1,605.00	\$1,249.00
	5 - 100ul	171-LH-729230	\$1,605.00	\$1,249.00
	30 - 300ul	171-LH-729240	\$1,605.00	\$1,249.00

Terms & Conditions: All prices exclude GST. Sale Prices valid on listed products only. Offer ends 31st July 2024.

Order today or call 1800 069 161 for more information.


MECHANICAL PIPETTES


The Mline Pipette family


WHERE PRECISION MEETS ACCURACY

Sartorius' mechanical pipette family - the Mline - offers excellent ergonomics, performance and safety in manual pipetting. It is particularly designed to maintain high accuracy and precision in repetitive and long-lasting pipetting. The Mline is an advanced and ergonomic choice for quality-conscious users who perform long and repetitive manual pipetting and is recommended by health and safety officers all over the world.

- ▶ Excellent ergonomics to protect from injuries with the exceptionally light pipetting and tip ejection forces
- ▶ High accuracy and precision, also in repetitive, long-lasting pipetting
- ▶ Minimised risk of contamination with Safe-Cone Filters and full autoclavability
- ▶ Increased safety with volume lock preventing accidental volume changes while pipetting

	Mline single channel pipette	Product Code	Original Price	SALE PRICE
SAVE \$100	0.1-3ul	171-725010	\$429.00	\$329.00
	0.5-10ul	171-725020	\$429.00	\$329.00
	2-20ul	171-725030	\$429.00	\$329.00
	10-100ul	171-725050	\$429.00	\$329.00
	20-200ul	171-725060	\$429.00	\$329.00
	100-1000ul	171-725070	\$429.00	\$329.00
	500-5000ul	171-725080	\$429.00	\$329.00
	1-10ml	171-725090	\$429.00	\$329.00

	Mline 8 channel pipette	Product Code	Original Price	SALE PRICE
SAVE \$227	0.5 - 10ul	171-725120	\$996.00	\$769.00
	5 - 100ul	171-725130	\$996.00	\$769.00
	30 - 300ul	171-725140	\$996.00	\$769.00

	Mline 12 channel pipette	Product Code	Original Price	SALE PRICE
SAVE \$277	0.5 - 10ul	171-725220	\$1,206.00	\$929.00
	5 - 100ul	171-725230	\$1,206.00	\$929.00
	30 - 300ul	171-725240	\$1,206.00	\$929.00

Terms & Conditions: All prices exclude GST. Sale Prices valid on listed products only. Offer ends 31st July 2024.

SARTORIUS


Pathtech™
First Choice In Science

MAXI VOLUME LIQUID HANDLING

Midi Plus Pipette Controller

Midi Plus is a lightweight cordless pipetting controller, which allows aspiration from bottles and tubes without the arm and hand elevations required with serological or volumetric pipettes. It fits all commonly used glass or plastic pipettes but can also be used with Sartorius 5 ml and 10 ml disposable tips.

- ▶ Precise fingertip control of variable speed ranges
- ▶ Fits all commonly used 1 - 100 ml glass or plastic pipettes
- ▶ Hydrophobic autoclavable filter prevents from over aspiration
- ▶ Convenient fold-out bench stand supports the unit and pipette when not in use

	Midi Plus (inc. AC Adaptor)	Product Code	Original Price	SALE PRICE
SAVE \$109	1 - 100ml	171-710931	\$608.00	\$499.00

The Prospenser and Prospenser Plus

The Prospenser and Prospenser Plus are high quality bottle-top dispensers with excellent chemical resistance. They are safe, convenient, and customizable, making them perfect for a wide range of dispensing needs.





- ▶ Extremely high chemical and thermal resistance due to ceramic piston, glass cylinder, and highly resistant fluoroplastic parts
 - ▶ The incremental volume setup ensures accurate and repeatable results
 - ▶ Fully autoclavable
 - ▶ Wide range of adapters included to fit most common bottle sizes
- The Prospenser Plus also includes 3 advanced features:**
- ▶ Media Recirculation System, Smart 360° Rotating Dispensing Head and Universal Connector for drying tubes and filters.

	Bottle Top Dispenser Prospenser	Product Code	Original Price	SALE PRICE
SAVE UP TO \$128	0.2 - 1ml	171-LH-723060	\$419.00	\$319.00
	0.4 - 2ml	171-LH-723061	\$419.00	\$319.00
	1 - 5ml	171-LH-723062	\$419.00	\$319.00
	2 - 10ml	171-LH-723063	\$419.00	\$319.00
	5 - 30ml	171-LH-723064	\$587.00	\$459.00
	10 - 60ml	171-LH-723065	\$608.00	\$499.00

	Bottle Top Dispenser Prospenser Plus	Product Code	Original Price	SALE PRICE
SAVE UP TO \$194	0.2 - 1ml	171-LH-723070	\$555.00	\$429.00
	0.4 - 2ml	171-LH-723071	\$555.00	\$429.00
	1 - 5ml	171-LH-723072	\$555.00	\$429.00
	2 - 10ml	171-LH-723073	\$555.00	\$429.00
	5 - 30ml	171-LH-723074	\$839.00	\$699.00
	10 - 60ml	171-LH-723075	\$923.00	\$729.00

Terms & Conditions: All prices exclude GST. Sale Prices valid on listed products only. Offer ends 31st July 2024.

Pipetting Accessories

 SAVE \$25	Linear stand for all Sartorius pipettes	Part # 171-725620	Original Price \$97.00	SALE PRICE \$72.00
 SAVE \$37	Carousel stand for 6 Mech pipettes	Part # 171-LH-725630	Original Price \$146.00	SALE PRICE \$109.00
 SAVE \$99	Electronic Single Pipette Charging stand + AC adaptor	Part # 171-730981	Original Price \$398.00	SALE PRICE \$299.00
 SAVE \$208	Electronic Pipette Charging carousel + AC adaptor	Part # 171-730991	Original Price \$807.00	SALE PRICE \$599.00

Terms & Conditions: All prices exclude GST. Sale Prices valid on listed products only. Offer ends 31st July 2024.

SARTORIUS

 **Pathtech**[™]
First Choice In Science

Pathtech provides calibration, maintenance, and repair services to extend the lifetime and ensure the continuous accuracy and precision of your pipettes via our Service Laboratories.

Why is Calibration and Maintenance Needed?

All pipettes manufacturers recommend that regular maintenance and calibration is applied to maintain reliable pipetting results and to maximize the lifetime of the pipette. Pipette calibration is a fundamental part of Good Laboratory Practice (GLP) and quality systems and must be considered a vital part of any laboratory regime where precise volumes of fluid need transferring or diluting.

Pipette performance is measured as accuracy and precision or how close the dispensed volumes are to the target and how close the results are to each other.

Why choose us for your next service?

- ▶ We service all major makes and models*
- ▶ We only use quality and reputable products
- ▶ Upfront costs, no hidden extras
- ▶ Experienced technicians
- ▶ Quick turnaround time
- ▶ Efficient customer service with the use of our Asset Management Portal
- ▶ 10 years plus accredited laboratory



Want to know more? Contact our friendly Customer Service Team today on 1800 069 161 or email info@pathtech.com.au for more information.

*We proudly calibrate and service most major models and makes of pipettes, such as Sartorius, Finnpiptette (ThermoFischer Scientific), Rainin/Mettler Toledo, Socorex, Axygen, Eppendorf, Gilson, Integra, Nichiryo (Lab Gear), Ovation and many more.

Call us now to place your order 1800 069 161
or visit pathtech.com.au for more info or to buy online

SARTORIUS

