

MyFuge® Mini Centrifuge with 2 Rotors

YOUR PERFECT LAB PARTNER



- › Clik-N-Lok™ INSTANT rotor attachment
- › 8 position microtube rotor & 4 position PCR strip rotor
- › Storage compartment for PCR rotor & 0.5ml adapters
- › Includes round & PCR rotors plus tube adapters
- › Near silent operation
- › Decelerates in just 1 second
- › Conserves Valuable Bench Space
- › Five colours

The MyFuge™ Mini centrifuge is a personal centrifuge that virtually fits in the palm of your hand. Ideal for quick spin downs of microtubes and PCR tubes, the MyFuge is extremely easy to use. Simply close the lid and the rotor instantly reaches 6,000 rpm. Open the lid, and the rotor quickly decelerates to a stop for removing samples.

An increased capacity, eight-position microtube rotor is included and comes installed. Also included is a second rotor for 0.2ml PCR tubes and strips which is stored (together with the rotor removal key) in a convenient compartment on the bottom of the centrifuge. When the rotor is not in use, storage in this compartment prevents loss or misplacement of the rotor.

For more information about the MyFuge Minis or any of our Centrifuge range please contact your Account Manager or our friendly Customer Service team on 1800 069 161 or via email info@pathtech.com.au

Ordering Information:

Product Code	Product Description
032-C1008-B	MyFuge™ Mini Centrifuge, Blue lid, with 2 rotors
032-C1008-C	MyFuge™ Mini Centrifuge, Clear lid, with 2 rotors
032-C1008-G	MyFuge™ Mini Centrifuge, Green lid, with 2 rotors
032-C1008-P	MyFuge™ Mini Centrifuge, Purple lid, with 2 rotors
032-C1008-R	MyFuge™ Mini Centrifuge, Red lid, with 2 rotors

Specifications:

Speed	6000 rpm / 2000 xg
Capacity	8 x 1.5/2.0 ml tubes 32 x 0.2 ml PCR tubes 4 x PCR Strips
Dimensions	14.5 x 17.8 x 11.4 cm
Weight	1kg
Electrical	Universal Voltage 100 to 240V; 50-60Hz, 20W
Warranty	2 Years



Storage Compartment



Includes Both Rotors