

Midi Plus

Pipette Controller 1-100ml

Easy and ergonomic pipetting of large volumes from bottles or tubes

Key Features

- › Ergonomic handle design with electronic operating buttons for reduced finger force
- › Stepless speed control for accurate dispensing
- › Possibility of dispensing with the force of gravity for slow, drop-by-drop dispensing (e.g., into culture media dish so as not to break the fine surface of the medium)
- › Fold-out stand saves bench space when Midi Plus is not in use and keeps the pipette detached from the table surface to avoid contamination
- › Works with all brands of Serological Pipettes 1-100ml
- › Low battery warning



Product Code: 171-710931

For more information about the Midi Plus or any of the Sartorius range please contact your Account Manager or our friendly Customer Service team on 1800 069 161 or via email info@pathtech.com.au

SARTORIUS

 **Pathtech**[™]
First Choice In Science